
Instruction for use

EN

1. Intended Purpose

DRYUM is a sterile single-use disposable medical device intended to deliver constant negative pressure to closed post-operative wounds, managing the microenvironment.

2. Intended User

DRYUM must be used by healthcare professionals who are properly trained, qualified, or certified in accordance with national legislation, following these instructions for use.

3. Indications

DRYUM is indicated for patients who benefit from negative pressure wound therapy, as it promotes healing by managing the wound environment and removing exudate and infectious fluids, thereby reducing the risk of infection and maintaining consistent moisture levels. Suitable wound types include closed surgical incisions.

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


b Soft connection
Connessione morbida


c Patch
Cerotto

2 Valve
Valvola

3 Adhesive strips x6
Strisce adesive x6

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h SAP
SAP

i Bottom cover
Cover

DRYUM is a single-use medical device that applies mechanically generated negative pressure to surgical incisions and post-operative wounds. It enhances wound healing by reducing post-operative complications and shortening healing time. Healing is promoted by removing exudate and infectious fluids and increasing blood flow to the wound site. **DRYUM** consists of a suction unit, two dressings, two valves, and six adhesive strips.

The dressing is made of a patch, a soft connector, and a connection tube. The patch consists of a hydrophobic fabric in direct contact with the wound and a skin-friendly silicone adhesive. The soft connector, containing an anti-collapse structure, allows fluid to flow between the patch and the connection tube. Its design minimizes pressure injury risk, avoids kinking of the tubing, and improves patient comfort. The connection tube links the soft connector to the suction unit through a valve. The suction unit generates negative pressure via a piston mechanism and can hold up to 55 ml of exudate, which is gelled by a superabsorbent polymer (SAP). The adhesive strips can be applied to the patch edges to improve adherence and prevent air leaks.

DRYUM is intended for a usage period of seven (7) days.

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5. Composition

Dressing: Polypropylene, Silicone, Polyurethane.
Suction Element: Polypropylene.
Adhesive strip: Polyurethane, acrylic adhesive.

6. Packaging and Sterility

The product comes in a box containing two pouches with dressings and one pouch with the suction unit. The product is sterile, sterilized via ethylene oxide. The product is NOT re-sterilizable.

DRYUM is available with the following dressing sizes:

12 cm x 20 cm
12 cm x 29 cm

7. Disposal

Dispose of the product according to the healthcare facility's waste disposal procedures and local regulations.

8. Contraindications

Negative pressure therapy must NOT be used in the following cases:

- Patients with known sensitivity or allergy to silicone or acrylic adhesives
- Malignant wounds (wound bed and/or wound edges) (except palliative care to improve quality of life)
- Wounds with confirmed and untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Necrotic wounds or eschar
- Wounds with exposed arteries, veins, nerves, or organs
- Anastomotic sites
- For emergency airway aspiration
- As a pleural, mediastinal, or thoracic drainage tube
- Burns, including partial-thickness burns

9. Warnings

- Hemostasis must be achieved before applying the dressing;
- Always ensure that the suction unit and tubing:
 - are not positioned in a way that could cause pressure injury to the patient;
 - are not dragged on the floor or pose a tripping or contamination hazard;
 - do not pose a strangulation risk or act as a tourniquet;
 - do not rest on or cross heat sources;
 - are not twisted or trapped under clothing or bandages, which may block negative pressure;
- The device contains small parts that could pose a choking hazard. Keep out of reach of children;

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- DRYUM dressings must only be applied and replaced by healthcare professionals. Patients should not remove the dressing;
- Each **DRYUM** dressing is intended for use on a single wound;
- The calibration marks on the suction unit are not clinical indicators. They inform the user of the correct functioning of **DRYUM** and/or the need to replace the dressing;
- Do not use DRYUM after the expiration date on the packaging;
- Do not use DRYUM if the packaging is damaged, as it may compromise sterility;
- Do not modify or replace any device components;
- Do not expose the product to direct sunlight;
- Store the product in a dry place;
- **DRYUM** is for single use only. Using any component on more than one patient may cause cross-contamination and infection;
- Do not use oily substances (ointments, creams) that may impair the patch's adhesion;
- Do not cut the patch, as this may compromise consistent negative pressure.

10. Precautions

- Keep the suction unit accessible for patient and clinician to verify correct function;
- **DRYUM** delivers negative pressure when the red piston is fully inside the suction unit. If the piston is visible, reload the device by pressing it again;
- Remove any hair from the application site that could compromise the adhesion of the dressing;
- Clothing may be worn after applying **DRYUM**. Ensure it does not interfere with tubing or block pressure flow;
- **DRYUM** may be used for up to seven (7) days.
- Always ensure that the dressing is positioned centrally with only the 3D spacer over and that the adhesive part is in contact with intact skin
- Do not cut the wound dressing as it may lead to loss of constant negative pressure

11. Adverse Event

No adverse events reported, as the product is not yet on the market.

Excessive bleeding is a serious risk associated with negative pressure wound therapy and may cause death or serious injury. Careful patient selection is critical, considering the contraindications, warnings, and precautions.

Monitor the wound and dressing for any signs of blood loss.

Monitor the patient for sudden changes in volume or color of exudate.

There is a potential risk of increased skin blistering, especially in patients with fragile or sensitive

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periwound skin. This may result from excessive shear forces during dressing application or removal, or prolonged tension at the skin-adhesive interface. The user or patient must report any adverse events related to the product to MEDICUD S.r.l at the following email: qualityassurance@medicud.it, and to the relevant national competent authority.

The product must only be applied, changed, and removed by trained healthcare professionals.

12. DRYUM application

1. Remove excess hair near the wound, if necessary disinfect and dry the area;
2. Open the cardboard box;
3. Take out the two pouches with dressings and store one in the patient file for the dressing change halfway through treatment;
4. Open the dressing pouch horizontally from the triangular end halfway to prevent valve or adhesive strips from falling out;
5. Remove the dressing without dropping the valve;
6. Remove the first liner (FIG.1) and place the dressing centered over the wound, ensuring only the central hydrophobic fabric is in contact with the wound, and stick the silicone adhesive to the skin, avoiding folds or wrinkles. The dressing should be positioned with the soft connector facing upwards if the wound is on the lower body, and facing downwards if on the upper body;
7. Remove the second and third liners (FIG.2) and ensure the patch adheres well to the skin (FIG.3). If misaligned or wrinkled, the dressing can be repositioned;
8. Apply the adhesive strips along the patch edges (FIG.4);
9. Once the dressing is in place, remove the tubing from its protective pouch. You may shorten the tubing with scissors based on the incision location;
10. Remove the suction unit pouch from the box;
11. Open the suction unit pouch;
12. Take the valve from the dressing pouch and screw it onto the suction unit (FIG. 5, FIG. 6, FIG. 7);
13. Insert the valve into the tube (FIG. 8) and ensure a proper seal;
14. Load the device in vertical position (red piston facing up) by pressing the piston with your thumb until it stops (FIG. 9 and FIG. 10), ensuring the tube is not kinked. Release the thumb and check that the piston remains in place. If it moves after a few seconds, repeat.

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13. DRYUM change

Change the dressing 3 to 4 days after initial application;

Leave the dressing connected to the suction unit, Reload the suction unit (only if air is in the reservoir) by pressing the red piston and checking the fill level using the

graduated scale. If the level is ≥ 40 ml (yellow), dispose of the entire DRYUM device (FIG.11). If the level is <40 ml (green), continue using the same suction unit and dispose only the dressing (FIG.12). Disconnect the valve from the suction unit and place it on a sterile surface. Removal and disposal of the dressing as biological waste.

NOTE: Once the valve is disconnected from the suction unit, do not press the red piston (FIG.13). The suction unit can only be reloaded when connected to a dressing via the check valve (FIG.14);

3. Take the new dressing from the patient file and repeat steps starting from point 4 in the "DRYUM Application" section.

14.DRYUM removal and disposal

1. Remove the second dressing 3/4 days after the first change;
2. Disconnect the valve from the suction unit, remove the dressing from the patient, and dispose of both the dressing and suction unit as biological waste.

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